
# InnerScan

## Body Composition Monitor

Model : BC-582

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## Eng Table of Contents

Introduction	1
Features and functions	1
Preparations before use	2
Getting accurate readings	2
Automatic shut-down function	2
Setting (resetting) time	3
Setting and storing personal data in memory	3
Switching the weight mode	4
What is the FITPLUS mode?	5
Inputting your menstrual period start day according to your natural cycle	6
Reminder to program your menstruation period start day	6
Taking weight and body fat readings	6
Programming the guest mode	7
Taking weight measurement only	7
Troubleshooting	8

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## INTRODUCTION

Thank you for selecting a Tanita Body Composition Monitor. This monitor is one in a wide range of home healthcare products produced by Tanita.

This Instruction Manual will guide you through the setup procedures and outline the button features. Additional information on healthy living can be found on our website [www.tanita.eu](http://www.tanita.eu). Tanita products incorporate the latest clinical research and technological innovations. All data is monitored by our Tanita Medical Advisory Board to ensure accuracy.

**Note:** Read this Instruction Manual carefully and keep it handy for future reference.

### ▲ Safety Precautions

Note that people with an electronic medical implant, such as a pacemaker, should not use a Body Composition Monitor as it passes a low-level electrical signal through the body, which may interfere with its operation. Pregnant women should use the weight function only. All other function are not intended for pregnant women.

Do not place this monitor on slippery surfaces such as a wet floor.

Please do not use scale immediately after bathing or strenuous exercise. Scale platform will be slippery when wet. In addition, droplets of water or sweat may activate the switches on the display. Please wipe platform with clean, dry cloth to remove any liquids.

### Important Notes for Users

This Body composition monitor is intended for adults and children (ages 5-17) standard and athletic body types.

**Tanita defines "athlete" as a person involved in intense physical activity of approximately 10 hours per week and who has a resting heart rate of approximately 60 beats per minute or less. Tanita's athlete definition includes "lifetime of fitness" individuals who have been fit for years but currently exercise less than 10 hours per week.**

**The body fat monitor function is not intended for pregnant women, professional athletes or bodybuilders.**


Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges. Tanita takes no responsibility for any kind of loss caused by the loss of recorded data.

Tanita takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.

**Note:** Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

## FEATURES AND FUNCTIONS



### Measuring Platform


1. FITPLUS Indicator
2. Arrow (Select) Buttons
3. Set Button
4. FITPLUS Button
5. Weight-Only/Off Button
6. Display Screen

### Accessories


AA-Size Batteries (4)  
(Included)

## PREPARATIONS BEFORE USE

Open as shown.

Direction of the batteries.

### Inserting the Batteries

Turn the scale over and remove the battery cover by releasing the latch, insert 4 x AA batteries (included) with the correct polarity as indicated.

**Note:** Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage. Please note that the included batteries from the factory may have decreased energy levels over time.

### Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

**Note:** To avoid possible injury, do not step on the edge of the platform.

## Handling Tips

This monitor is a precision instrument utilising state-of-the-art technology. To keep the unit in the best condition, please follow these instructions carefully:

- Do not attempt to disassemble the measuring platform.
- Store the unit horizontally, and place it so that the buttons will not be pressed accidentally.
- Avoid excessive impact or vibration to the unit.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner (applied to a cloth first) to keep them shiny; avoid using soaps.
- Do not step on the platform with wet feet.
- Do not drop any objects onto the platform.

## GETTING ACCURATE READINGS

Heels centered on electrodes

Toes may overhang measuring platform

**Note:** An accurate reading will not be possible if the soles of your feet are not clean, or if your knees are bent or you are in a sitting position.

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping onto the measuring platform. Be sure your heels are correctly aligned with the electrodes on the measuring platform. Don't worry if your feet appear too large for the unit, accurate readings can still be obtained if your toes overhang the platform. It is best to take readings at the same time of day. Try to wait about three hours after rising, eating, or exercising before taking measurements. While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage of change as long as the readings are taken in a consistent manner. To monitor progress, compare weight and body fat percentage taken under the same conditions over a period of time.

## AUTOMATIC SHUT-DOWN FUNCTION

The automatic shut-down function shuts off the power automatically in the following cases:


- If you interrupt the measurement process. The power will shut down automatically within 10 to 20 seconds, depending upon the type of operation.
- If an extreme weight is applied to the platform.
- During programming, if you do not touch any of the buttons within 60 seconds.\*
- After you have completed the measuring process.

\***Note:** If the power shuts off automatically, repeat the steps from "Setting and Storing Personal Data in Memory" (p.3).

## SETTING (RESETTING) TIME


### 1. Set Hour

Press and hold for 3 seconds to turn on the unit. Press to set hour and then press . The unit will beep once to confirm.

### 2. Set Minute

Press to set minute and then press . The unit will beep twice to confirm the programming. The unit will then shut down automatically.

### Resetting time

Press and hold for 3 seconds to turn on the unit. Then follow the procedure 1 and 2.

**Note:** If you chose FitPlus mode, please re-input your menstrual period start date. Select your personal number and enter how many days ago your menstrual period started.


## SETTING AND STORING PERSONAL DATA IN MEMORY

You will receive readings only if data has been programmed into one of the personal data memories or the Guest mode. The monitor can store personal details for up to 5 people.

### 1. Turn on the Power


Press to turn on the unit. The unit will beep to confirm activation, the Personal Button numbers will be displayed, and the display will flash.

**Note:** If you make a mistake or want to turn the unit off before you have finished programming it, press .


### 2. Select a Personal Key Number

Press to select a Personal Button. Once you reach the Personal Button number you wish to use, press . The unit will beep once to confirm.

### 3. Set Age

Use to scroll through the numbers. When you reach your age, press . The unit will beep once to confirm.

Continue to the next page.


#### 4. Select Female or Male

Press to scroll through Female () or Male () , Female/Athlete () / ) and Male/Athlete () / ) settings, then press . The unit will beep once to confirm.

#### 5. Specify the Height

Use to specify Height and then press . The unit will beep once to confirm.

#### 6. Select Activity Level

To use the DCI (Daily Calorie Intake) feature properly, you need to set the Activity Level. Use to specify Activity Level. Press to confirm.

- Level 1 = Inactive lifestyles (little or no exercise)
  - Level 2 = Moderately active lifestyles (occasional, low intensity exercise)
  - Level 3 = Adult involved in intense physical activity (see p.1. for Janita's definition of Athlete)
- Note: If you select the Athlete Mode, it will automatically set the Activity Level to Level 3.**


The unit will beep twice and the display will flash all data (Male/Female, Age, Height, Activity Level, etc.) three times to confirm the programming. The power will then shut down automatically.

#### 7. Programming FITPLUS Mode (Women Only)

After selecting Activity Level, use to select "on" or "OFF" of the FITPLUS Mode, and press . The unit will beep once to confirm. The indication () appears on the upper right side of the display to indicate that you have selected the FITPLUS Mode.

After selecting "on", use to enter how many days ago your menstrual period started and press to confirm. The unit will beep once to confirm. The unit will then beep twice and the display will flash all of the programmed user data three times to confirm the programming. The power will then shut down automatically.

### SWITCHING THE WEIGHT MODE


Select your preferred weight mode using the weight mode switch on the bottom of the scale.

### WHAT IS THE FITPLUS MODE?

Women have an optimal fitness phase, ovulation window and days where their body naturally needs to rest. FITPLUS mode gives you healthy prompts according to your natural cycle, these are indicated using an easy-to-read traffic light system.

- 'GREEN' days, around days 5 to 14 in the Follicular Phase is when exercising is optimal with around double the calories being burnt than on other days
- 'YELLOW' days, around days 13 to 16 indicate likely ovulation
- 'RED' days, from days 17 to 28 in the Luteal Phase should be less active as the body naturally needs more rest

Measurement Results in FITPLUS Mode

Indicator color

- Green= Good diet period
- Yellow= Changing period
- Red= Relax period

#### Measurement Results in FITPLUS Mode

When you have programmed the FITPLUS Mode, the FITPLUS Mode lamp is illuminated while all measurements are displayed.

**Note:** The light does not turn on in the cases below:

- Have no menstrual period
- Menopausal
- Pregnant
- Under 18 years old
- No programming of menstrual period start day

## INPUTTING YOUR MENSTRUAL PERIOD START DAY ACCORDING TO YOUR NATURAL CYCLE



To ensure you receive personalised readings your menstrual period start day needs to be input every month.

### Program your menstrual period start day.

Press to turn on. Use to select your Personal data number. Press to confirm.

The unit will beep and the display will show the programmed data.

The unit will beep again and the display will show "0.0". Press . Use to set the menstrual period start day. Press to confirm. The power will then shut down automatically.



## REMINDER TO PROGRAM YOUR MENSTRUATION PERIOD START DAY

The lamp flashes to remind you to set the first day of your menstrual cycle.

Program the menstrual period start day while the lamp is flashing: The FITPLUS Mode lamp will not turn on again until the programming is executed.

## TAKING WEIGHT AND BODY FAT READINGS


### 1. Select Personal Data Number

Press to turn on. Use to select your Personal data number. Press to confirm. The unit will beep and the display will show the programmed data.

**The unit will beep again and the display will show "0.0".  
Now step onto the platform.**

**Note:** If you step onto the platform before "0.0" appears the display will show "Error" and you will not obtain a reading. Furthermore, if you do not step onto the measuring platform within 30 seconds after "0.0" appears, the power will shut-off automatically.

### 2. Obtain Your Readings


Your weight will be shown first. Continue to stand on the platform. "00000" will appear on the display and disappear one by one from left to right. Your Body Fat percentage (%FAT) and Body Fat Range will appear on the display screen.

**Note: Do not step off until Body Fat % is shown.**

**Note: For children (age 5-17), the unit will only display the Weight and the Body Fat %.**

To see other readings, **carefully step off scale**, press continually to display readings in the order listed below. (Press for reverse.)

will not function while remaining on scale.

### 3. Recall Function

The previous readings of the following features can be obtained using the recall function.



- Weight
- Body Fat %
- Total Body Water %
- Muscle Mass/Physique Rating
- DCI/Metabolic Age
- Bone Mass
- Visceral Fat Level

To obtain previous readings, press while the current readings are being displayed.

The display will show the previous readings. To see other previous readings, press for each desired reading. Press again to go back to the current result display.

**Note:** For Children (age 5-17), the unit will only display the previous Weight and the Body Fat %.

## PROGRAMMING THE GUEST MODE


Guest mode allows you to programme the unit for a one-time use without resetting a Personal data number.

Press to activate the unit and select 'G'. Then follow the steps: set age, select female or male, and specify height. Once programmed please press set to confirm. The unit will beep, then flash through the data and go to "0.0".

**Note:** If you step onto the platform before "0.0" appears the display will show "Error" and you will not obtain a reading.

## TAKING WEIGHT MEASUREMENT ONLY

### Select Weight-Only Key

Press . After 2 or 3 seconds, "0.0" will appear in the display and the unit will beep. When "0.0" appears in the display, step onto the platform. The body weight value will appear. The power will then shut down automatically.

## TROUBLESHOOTING

### If the following problems occur ... then...

- **[Lo]** appears on the display, or all the data appears and immediately disappears.  
Batteries are low. When this message appears, be sure to replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new batteries.  
Note: Your settings will not be erased from the memory when you remove the batteries.
- **[Error]** appears while measuring.  
Please stand on the platform keeping movement to a minimum.  
The unit cannot accurately measure your weight if it detects movement.  
The Body Fat Percentage is more than 75%, readings cannot be obtained from the unit.
- **The Body Fat Percentage reading does not appear or [.....] appears after the weight is measured.**  
Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned with the guides on the measuring platform.
- **[UL]** appears while measuring.  
Reading cannot be obtained if the weight capacity is exceeded.

### Specifications

BC-582

### Weight Capacity

330lb (150kg) (23st 8lb)

### Weight Increments

0.2lb (0.1kg) (0.2lb)

### Body Fat Increments

0.1%

### Personal Data

5 memories

### Power Supply

4 AA batteries (included)

### Power Consumption

60 mA maximum

### Measuring Current

50 kHz, 6.25 kHz, 100 µA

The device features radio interference suppression in compliance with EC regulation 2004/108/EC

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